DIGITAL SOFTWARE NEWS

FOR THE

PDP-8 & PDP-12

Software Information Service Digital Equipment Corporation Maynard, Massachusetts 01754

DIGITAL SOFTWARE NEWS FOR THE PDP-8 AND THE PDP-12

Digital Software News for the PDP-8 & PDP-12 is designed to provide users of PDP-8 and PDP-12 software with up to date information on software problems, programming notes, and new and revised software. Each article, other than announcements and documentation corrections, is coded sequentially by system program in the lower right corner for ease in maintaining a quick reference file. The original material for such a file is supplied in the Software Performance Summary for the PDP-8 This publication which is placed in each basic software kit, is a collection of all current information on known software problems, patches, and programming notes. new versions correct software problems and obsolete patches, and reprinted manuals include programming notes and manual corrections, new articles in the Digital Software News will announce the revised software and specify by code which article should be removed from your Software Performance Summary file. Articles may also be replaced when new information becomes available; such as, procedure to circumvent a problem may replace the original report of the problem. All articles in the Digital Software News which should be added to your file will be appropriately coded.

Any questions or problems on the articles contained in these publications or concerning the use of Digital's software should be reported to the Software Specialist or Sales Engineer at your nearest Digital Office.

Additional copies of the current <u>Software Performance Summary</u> and updated issues of the PDP-8 and PDP-12 software price lists are available at no charge upon request from the Program Library. As with usual orders for software, the Program Library will accept orders received directly from U.S. customers. All other customer orders must be routed through the nearest Digital Representative.

	ø
	٠

PDP-8/12 DIGITAL SOFTWARE NEWS

OCTOBER-NOVEMBER

TABLE OF CONTENTS

SOFTWARE PROBLEMS

PDP-8 AIPOS MONITOR Problem on RUBOUT	AIPOS MONITOR	1
FOCAL 1969 (DEC-Ø8-AJAE Optimum use of the PDP-8/E	FOCAL	17
PDP-12 MASS SPECTROSCOPY HANDLER (MASH) Limitations in Mash	MASH	2

		4) 4
		\$
		6

AIPOS MONITOR

Problem on RUBOUT

In the AIPOS Monitor, on RUBOUT, the command interpreter sometimes stores information in the wrong file descriptor block. To correct this problem, the following binary patch should be implemented.

Block 1 (add 1 to the starting block) of JOB CONTROL

LOCATION	OLD CONTENTS	NEW CONTENTS
224	7 Ø Ø 1	7 Ø Ø Ø

To implement a source patch, $CM\emptyset2$, line 417 [IAC] should be deleted. The segment of the source should now read:

TAD I P1FAUX AND P7 IAC DCA FILCNT JMP CMI18Ø

			100
			ž
			ŵ
			ž.

FOCAL 1969 (DEC- \emptyset 8-AJAE)

Optimum use of the PDP-8/E

The following patch uses the KCF instruction of the PDP-8/E to eliminate input buffer overflow. The patch should be assembled with PAL III so that the binary can be loaded easily each time FOCAL is loaded.

```
/PATCH TO FOCAL 1969 (AJAE)
               /TO CURE BUFFER OVERFLOW PROBLEM (PDP-8/E ONLY)
               /MAKES USE OF EXTRA 8/E KEYBOARD IOT (KCF)
               /SO MAKING IT IMPOSSIBLE FOR OVERFLOW TO OCCUR
               /PATCHES FOCAL INTERRUPT PROCESSOR
               /LISTING PAGE 44
               /AND FOCAL TTY INPUT ROUTINE
               /LISTING PAGE 45
                OPATCH=2414
                INBUF = 34
                P177 = 106
                MBREAK= 2602
                C200 = 123
                EXIT= 2646
                RECOVR= 2740
                SIN = 2662
                INRET = 2675
                KCF = 6030
                *2627
               /PATCH TO INTERRUPT PROCESSOR
                                          /INPUT CHARACTER
2627
       6034
                       KRS
                                          /CLEAR FLAG
                       KCF
2630
       6030
                                          /IGNORE BIT 8
                       AND P177
2631
       0106
                                          /BLANK?
       7450
                       SNA
2632
                                          /YES-GO INITIATE NEXT READ
2633
       5245
                       JMP EXIT-1
                                          /SET BIT 8 ON
2634
       1123
                       TAD C200
                       DCA SIN
2635
       3262
                       TAD SIN
       1262
2636
                       TAD MBREAK
2637
       1202
                                          /CTRL C?
                       SNA CLA
2640
       7650
2641
       5340
                       JMP RECOVR
                                          /YES
                       TAD SIN
2642
       1262
                                          /STORE CHARACTER
                       DCA INBUF
2643
       3034
                        SKP
2644
       7410
                                          /INITIATE NEXT READ- CHAR. WAS A BLANK
2645
       6032
                       KCC
               *2673
                        JMP I .+1
                                          /PATCH TO INPUT ROUTINE
2673
       5674
                        OPATCH
2674
       2414
               *OPATCH
                                          /CLEAR INPUT BUFFER
2414
       3034
               OPATCH
                       DCA INBUF
                                          /INITIATE NEXT READ
                        KCC
2415
       6032
                        TAD I XOUTL
                                          /GET CHARACTER
2416
       1621
                                          /RETURN TO INPUT ROUTINE
2417
       5620
                        JMP I .+1
                        INRET
2420
       2675
                        2676
2421
       2676
               XOUTL,
               $
```

PDP-12 MASS SPECTROSCOPY HANDLER (MASH)

Limitations in MASH

The current library version of MASH assumes a non-quadrupole mass spectrometer and a 12-bit A/D as part of its standard hardware configuration. The following patches eliminate those restrictions.

A. For a 15-bit A/D:

BLOCK	LOCATION	NEW CONTENTS	OLD CONTENTS	SOURCE
33	122	3 2 2 Ø	3221	/RESCALE PEAK
	123	1214	1215	/X-Y PARA.
	134	3220	3221	

B. For a quadrupole mass spectrometer of the Finnegan variety:

BLOCK	LOCATION	NEW CONTENTS	OLD CONTENTS	SOURCE
6	217	15 Ø 1	22Ø5	/RESOLUTION IS
	22Ø	2323	2317	/NOW MASS RANGE
	221	4 Ø 2 2	1425	
	222	Ø116	2411	
	223	Ø7Ø5	1716	
	315	2225	22Ø1	
21	331	2	7776	/FINNEGAN PEAK
	332	33ØØ	2365	/WIDTH CONSTANT
	333	Ø	6 Ø 5 Ø	
2 5	15Ø	5	4	/FINNEGAN RELAY
	1 51	4	5	/PROBLEMS
26	166	Ø	-	/WAITM, Ø/WAITING /RELAY 2
	167	73ØØ	-	/CLA CLL
	17Ø	3Ø15		/DCA 15
	171	115Ø	-	/TAD STSCN1
	172	4563	***	/JMS I DOSCN1
	173	2015	-	/ISZ 15
	174	5371	-	/JMP WAITM+3
	175	5766	-	/JMP I WAITM
2 7	46	1376	1150	/LOAD RELAY 2
	176	2		
30	13	4775	7 3 Ø Ø	/JMS WAITM
2.4	175	566	8000	A TITOMII CONOMANT
34	354	0	1	/PEAK WIDTH CONSTANT
52	362	4 Ø	1031	/ELIMINATE 50 UNIT
	363	4 Ø	1752	/MASS CHECK

To Replace